

Pragmin Antibody

Catalog No: #45379

Package Size: #45379-1 50ul #45379-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pragmin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Pragmin Antibody detects endogenous levels of total Pragmin
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Pragmin
Target Name	Pragmin
Other Names	SGK223, Homolog of rat pragma of Rnd2, PRAGMIN, Sugen kinase 223, Tyrosine-protein kinase SgK223
Accession No.	Swiss-Prot#: Q86YV5NCBI Gene ID:157285
Uniprot	Q86YV5
GeneID	157285;
Calculated MW	149 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only