

KIAA0494 antibody

Catalog No: #23101



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	KIAA0494 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 243 and 481 of KIAA0494
Target Name	KIAA0494
Accession No.	Swiss-Prot:O75071Gene ID:9813
Uniprot	O75071
GeneID	9813;
Concentration	0.6mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 55kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250
Immunofluorescence: 1:100-1:200


Images




Sample (30 ug of whole cell lysate)
A: Molt-4
B: Raji
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded A549 xenograft, using KIAA0494 antibody at 1: 100 dilution.



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using KIAA0494 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only