

Tim17 antibody

Catalog No: #23106

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Tim17 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 108 and 121 (Q99595) of TIM17
Target Name	Tim17
Other Names	TIM17; TIM17A
Accession No.	Swiss-Prot:Q99595Gene ID:10440
Uniprot	Q99595
GeneID	10440;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 18kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: Hela
12% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only