

OR5H6 Antibody

Catalog No: #45447

Package Size: #45447-1 50ul #45447-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR5H6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OR5H6 Antibody detects endogenous levels of total OR5H6
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR5H6
Target Name	OR5H6
Other Names	Olfactory receptor 5H6; Olfactory receptor OR3-11; OR5H6
Accession No.	Swiss-Prot#: Q8NGV6NCBI Gene ID:79295
Uniprot	Q8NGV6
GeneID	79295;
Calculated MW	34 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis OR5H6 using Jurkat whole cell lysates

Note: This product is for in vitro research use only