

ANP32D Antibody

Catalog No: #45479

Package Size: #45479-1 50ul #45479-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ANP32D Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	ANP32D Antibody detects endogenous levels of total ANP32D
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ANP32D
Target Name	ANP32D
Other Names	PP32R2
Accession No.	Swiss-Prot#: O95626NCBI Gene ID:23519
Uniprot	O95626
GeneID	23519;
Calculated MW	15 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis ANP32D using COLO205 whole cell lysates

Note: This product is for in vitro research use only