

KR1_HHV11 Antibody

Catalog No: #45484

Package Size: #45484-1 50ul #45484-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KR1_HHV11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Human herpesvirus 1
Specificity	KR1_HHV11 Antibody detects endogenous levels of total KR1_HHV11
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from KR1_HHV11
Target Name	KR1_HHV11
Other Names	Serine/threonine-protein kinase US3
Accession No.	Swiss-Prot#: P04413NCBI Gene ID:2703401
Uniprot	P04413
GeneID	2703401;
Calculated MW	53kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-2000

Note: This product is for in vitro research use only