

CT47A Antibody

Catalog No: #45530

Package Size: #45530-1 50ul #45530-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CT47A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CT47A Antibody detects endogenous levels of total CT47A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CT47A
Target Name	CT47A
Other Names	Cancer/testis antigen 47A Cancer/testis antigen 47 CT47
Accession No.	Swiss-Prot#: Q5JQC4NCBI Gene ID:100507170/255313/653282/728036/728042/728049/728062/728072/728075/728082/728090/728096
Uniprot	Q5JQC4
GeneID	100507170;255313;653282;728036;728042;728049;728062;728072;728075;728082;728090;728096;
Calculated MW	30 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis CT47A using HeLa whole cell lysates

Note: This product is for in vitro research use only