

FAM18B antibody

Catalog No: #23127



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM18B antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 58 and 119 of FAM18B
Target Name	FAM18B
Accession No.	Swiss-Prot:Q9NYZ1Gene ID:51030
Uniprot	Q9NYZ1
GeneID	51030;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 24kd
Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: HeLa
12% SDS PAGE
FAM18B antibody diluted at 1: 1000

Note: This product is for in vitro research use only