

ZNF225 Conjugated Antibody

Catalog No: #C47864



Package Size: #C47864-AF350 100ul #C47864-AF405 100ul #C47864-AF488 100ul
 #C47864-AF555 100ul #C47864-AF594 100ul #C47864-AF647 100ul
 #C47864-AF680 100ul #C47864-AF750 100ul #C47864-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	ZNF225 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Recognizes endogenous levels of ZNF225 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human ZNF225.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Zinc finger protein 225
Accession No.	Swiss-Prot#:Q9UK10NCBI Gene ID:7768
Calculated MW	82
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.