

TMED9 antibody

Catalog No: #23133

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TMED9 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 200 of human TMED9
Target Name	TMED9
Accession No.	Swiss-Prot:Q9BVK6Gene ID:54732
Uniprot	Q9BVK6
GeneID	54732;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 27kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images




Sample (30 ug of whole cell lysate)

A: H1299

12% SDS PAGE

Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded DLD-1 xenograft, using TMED9 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only