

PRPF39 antibody

Catalog No: #23136

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | PRPF39 antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Purified by antigen-affinity chromatography. |
| Applications | WB IHC |
| Species Reactivity | Hu |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 146 and 385 of PRPF39 |
| Target Name | PRPF39 |
| Accession No. | Swiss-Prot:Q86UA1Gene ID:55015 |
| Uniprot | Q86UA1 |
| GeneID | 55015; |
| Concentration | 0.7mg/ml |
| Formulation | Supplied in 1XPBS, 1%BSA, 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |


Application Details

Predicted MW: 78kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

Images




Sample (30 ug of whole cell lysate)

A: Raji

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using PRPF39 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only