

ABHD3 Antibody

Catalog No: #45742

Package Size: #45742-1 50ul #45742-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABHD3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ABHD3 Antibody detects endogenous levels of total ABHD3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ABHD3
Target Name	ABHD3
Other Names	ABHD3, LABH3, Lung alpha/beta hydrolase 3
Accession No.	Swiss-Prot#: Q8WU67NCBI Gene ID:171586
Uniprot	Q8WU67
GeneID	171586;
Calculated MW	46 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only