

RNA binding motif protein 28 antibody

Catalog No: #23138



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	RNA binding motif protein 28 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 475 and 715 of Human RBM28
Target Name	RNA binding motif protein 28
Accession No.	Swiss-Prot:Q9NW13Gene ID:55131
Uniprot	Q9NW13
GeneID	55131;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 86kd  
Western blotting: 1:500-1:3000  
Immunofluorescence: 1:100-1:200

Images



Sample (30 ug of whole cell lysate)  
A: 293T  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Confocal immunofluorescence analysis (Olympus FV10i) of paraformaldehyde-fixed HeLa, using RBM28 antibody (Green) at 1: 100 dilution and alpha-tubulin antibody (Red) at 1: 2500.

---

Note: This product is for in vitro research use only