

SLC4A2 Antibody

Catalog No: #45768

Package Size: #45768-1 50ul #45768-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC4A2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB IHC (peptide)
Species Reactivity	Hu Ms Rt
Specificity	SLC4A2 Antibody detects endogenous levels of total SLC4A2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SLC4A2
Target Name	SLC4A2
Other Names	SLC4A2, Anion exchange protein 2, AE 2, Anion exchanger 2, Anion exchanger 2 type a, Anion exchanger 2 type b1, Anion exchanger 2 type b2, BND3L, EPB3L1, AE2, HKB3, MPB3L, Anion exchanger type 2, Erythrocyte band 3, NBND3
Accession No.	Swiss-Prot#: P04920NCBI Gene ID:6522
Uniprot	P04920
GeneID	6522;
Calculated MW	137 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

IHC 1:200

Note: This product is for in vitro research use only