

ABCA3 Antibody

Catalog No: #45789

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCA3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ABCA3 Antibody detects endogenous levels of total ABCA3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ABCA3
Target Name	ABCA3
Other Names	ABC3
Accession No.	Swiss-Prot#: Q99758NCBI Gene ID:21
Uniprot	Q99758
GeneID	21;
Calculated MW	191 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only