

RAP2C antibody

Catalog No: #23146



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	RAP2C antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 161 of Human RAP2C
Target Name	RAP2C
Accession No.	Swiss-Prot:Q9Y3L5Gene ID:57826
Uniprot	Q9Y3L5
GeneID	57826;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 21kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250

Images



Sample(30 ug whole cell lysate)  
A: MOLT4  
12% SDS PAGE  
Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded N87 xenograft, using RAP2C antibody at 1: 100 dilution.

---

Note: This product is for in vitro research use only