

CPZ Antibody

Catalog No: #45852

Package Size: #45852-1 50ul #45852-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CPZ Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Rt
Specificity	CPZ Antibody detects endogenous levels of total CPZ
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CPZ
Target Name	CPZ
Other Names	CPZ, Carboxypeptidase Z, Metallo-carboxypeptidase Z
Accession No.	Swiss-Prot#: Q66K79NCBI Gene ID:8532
Uniprot	Q66K79
GeneID	8532;
Calculated MW	74 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only