

TAL1/2 (acetyl-Lys221+222/Lys36+37) Conjugated Antibody

Catalog No: #CHW211

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #CHW211-AF350 100ul #CHW211-AF405 100ul #CHW211-AF488 100ul

#CHW211-AF555 100ul #CHW211-AF594 100ul #CHW211-AF647 100ul

#CHW211-AF680 100ul #CHW211-AF750 100ul #CHW211-Biotin 100ul

Description

Product Name	TAL1/2 (acetyl-Lys221+222/Lys36+37) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu,Ms
Specificity	Recognizes endogenous levels of TAL1/2 (AcLys221+222/Lys36+37) protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human TAL1/2.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TAL1; BHLHA17; SCL; TCL5; T-cell acute lymphocytic leukemia protein 1; TAL-1; Class A basic helix-loop-helix protein 17; bHLHa17; Stem cell protein; T-cell leukemia/lymphoma protein 5; TAL2; BHLHA19; T-cell acute lymphocytic leukemia protein 2; TAL-2; Class A basic helix-loop-helix protein 19; bHLHa19
Accession No.	Swiss-Prot#:P17542; Q16559NCBI Gene ID:6886; 6887
Uniprot	P17542
GeneID	6886;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	34
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only