

## Preclinical Compound Library

Catalog No: #L3410

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Preclinical Compound Library
Brief Description	Preclinical Compound Library is a collection of 212 compounds that are in preclinical phase with clear targets and detailed information on disease indication and reference.
Storage	Powder or pre-dissolved DMSO solutions in 96 well plate with optional 2D barcodeShipped with blue ice; Stable for One year as powder, 6 months at - 20 ° C in DMSO or 12months at -80 ° C in DMSO

## Application Details

Number of Compounds:212

## Product Description

A unique collection of 212 preclinical compounds for high throughput screening (HTS) and high content screening (HCS); Hot lead compounds for new drug discovery and development; Detailed compound information with structure, target, activity, IC50 value, and biological activity description; Structurally diverse, medicinally active, and cell permeable; NMR and HPLC validated to ensure high purity and quality;

Note: This product is for in vitro research use only