

CDR1 Antibody

Catalog No: #45874



Package Size: #45874-1 50ul #45874-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CDR1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CDR1 Antibody detects endogenous levels of total CDR1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CDR1
Target Name	CDR1
Other Names	Cerebellar degeneration-related antigen 1 CDR34
Accession No.	Swiss-Prot#: P51861NCBI Gene ID:1038
Uniprot	P51861
GeneID	1038;
Calculated MW	31 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only