

C1orf165 antibody

Catalog No: #23150

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C1orf165 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 225 of C1orf165
Target Name	C1orf165
Accession No.	Swiss-Prot:Q7L4P6Gene ID:79656
Uniprot	Q7L4P6
GeneID	79656;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 48kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images




Sample (30 ug of whole cell lysate)


A: Molt-4

10 % SDS PAGE

Primary antibody diluted at 1: 1000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using C1orf165 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only