

CKLF7 Antibody

Catalog No: #45887

Package Size: #45887-1 50ul #45887-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CKLF7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CKLF7 Antibody detects endogenous levels of total CKLF7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CKLF7
Target Name	CKLF7
Other Names	CMTM7, CKLFSF7
Accession No.	Swiss-Prot#: Q96FZ5NCBI Gene ID:112616
Calculated MW	19 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis CKLF7 using K562 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.