

KIAA0101 Polyclonal Antibody

Catalog No: #27394



Package Size: #27394-1 50ul #27394-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	KIAA0101 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human KIAA0101 (NP_055551.1).
Other Names	PCLAF; KIAA0101; L5; NS5ATP9; OEATC; OEATC-1; OEATC1; PAF; PAF15; p15(PAF); p15/PAF; p15PAF; PCNA-associated factor
Accession No.	Swiss-Prot#:Q15004NCBI Gene ID:9768
Uniprot	Q15004
GeneID	9768;
Calculated MW	12kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using KIAA0101 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only