

ZSCAN4 Polyclonal Antibody

Catalog No: #27589



Package Size: #27589-1 50ul #27589-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZSCAN4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZSCAN4 (NP_689890.1).
Other Names	ZSCAN4; ZNF494; zinc finger and SCAN domain containing 4
Accession No.	Swiss-Prot#:Q8NAM6NCBI Gene ID:201516
Uniprot	Q8NAM6
GeneID	201516;
Calculated MW	49kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZSCAN4 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only