

LRSAM1 Antibody

Catalog No: #45954

Package Size: #45954-1 50ul #45954-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRSAM1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LRSAM1 Antibody detects endogenous levels of total LRSAM1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LRSAM1
Target Name	LRSAM1
Other Names	LRSAM1, CMT2P, HTAL, Tsg101-associated ligase, RIFLE, TAL
Accession No.	Swiss-Prot#: Q6UWE0NCBI Gene ID:90678
Uniprot	Q6UWE0
GeneID	90678;
Calculated MW	84 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis LRSAM1 using HT29 whole cell lysates

Note: This product is for in vitro research use only