

BTN3A3 Polyclonal Antibody

Catalog No: #28475



Package Size: #28475-1 50ul #28475-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	BTN3A3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human BTN3A3 (NP_008925.1).
Other Names	BTN3A3; BTF3; BTN3.3; butyrophilin subfamily 3 member A3
Accession No.	Swiss-Prot#:O00478NCBI Gene ID:10384
Uniprot	O00478
GeneID	10384;
Calculated MW	65kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using BTN3A3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only