

ENOX1 Polyclonal Antibody

Catalog No: #28480



Package Size: #28480-1 50ul #28480-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ENOX1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ENOX1 (NP_060463.2).
Other Names	ENOX1; CNOX; PIG38; bA64J21.1; cCNOX; ecto-NOX disulfide-thiol exchanger 1
Accession No.	Swiss-Prot#:Q8TC92NCBI Gene ID:55068
Uniprot	Q8TC92
GeneID	55068;
Calculated MW	73kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ENOX1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only