

PARP12 Polyclonal Antibody

Catalog No: #28501



Package Size: #28501-1 50ul #28501-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PARP12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PARP12 (NP_073587.1).
Other Names	PARP12; ARTD12; MST109; MSTP109; ZC3H1; ZC3HDC1; poly [ADP-ribose] polymerase 12
Accession No.	Swiss-Prot#:Q9H0J9NCBI Gene ID:64761
Uniprot	Q9H0J9
GeneID	64761;
Calculated MW	79kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PARP12 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only