

EIF3A Antibody

Catalog No: #45975

Package Size: #45975-1 50ul #45975-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EIF3A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	EIF3A Antibody detects endogenous levels of total EIF3A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human EIF3A
Target Name	EIF3A
Other Names	EIF3A, Centrosomin homolog, Cytoplasmic protein p167, EIF3S10, EIF3 p167, EIF3 p180, EIF-3-theta, EIF3 p185, EIF3, p180 subunit, KIAA0139, p167, p180, p185, TIF32, EIF3, EIF3-p170, EIF3-theta
Accession No.	Swiss-Prot#: Q14152NCBI Gene ID:8661
Uniprot	Q14152
GeneID	8661;
Calculated MW	167kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis EIF3A using RAW264.7 whole cell lysates



Western blotting analysis using EIF3A Antibody #45975.

Note: This product is for in vitro research use only