

RPUSD2 Polyclonal Antibody

Catalog No: #28550



Package Size: #28550-1 50ul #28550-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RPUSD2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RPUSD2 (NP_689473.1).
Other Names	RPUSD2; C15orf19; C18B11; RNA pseudouridylate synthase domain containing 2
Accession No.	Swiss-Prot#:Q8IZ73NCBI Gene ID:27079
Uniprot	Q8IZ73
GeneID	27079;
Calculated MW	Refer to Figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only