

## COA3 Polyclonal Antibody

Catalog No: #28632

Package Size: #28632-1 50ul #28632-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	COA3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human COA3 (NP_001035521.1).
Other Names	COA3; CCDC56; COX25; HSPC009; MITRAC12; cytochrome c oxidase assembly factor 3
Accession No.	Swiss-Prot#:Q9Y2R0NCBI Gene ID:28958
Uniprot	Q9Y2R0
GeneID	28958;
Calculated MW	12kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of MCF7 cells, using COA3 at 1:1000 dilution.

## Background

This gene encodes a member of the cytochrome c oxidase assembly factor family. Studies of a related gene in fly suggest that the encoded protein is localized to mitochondria and is essential for cytochrome c oxidase function.

---

Note: This product is for in vitro research use only