

ZNF296 Polyclonal Antibody

Catalog No: #28657



Package Size: #28657-1 50ul #28657-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF296 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZNF296 (NP_660331.1).
Other Names	ZNF296; ZFP296; ZNF342; zinc finger protein 296
Accession No.	Swiss-Prot#:Q8WUU4NCBI Gene ID:162979
Uniprot	Q8WUU4
GeneID	162979;
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF296 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only