

ZNF177 Polyclonal Antibody

Catalog No: #28709



Package Size: #28709-1 50ul #28709-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF177 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ZNF177 (NP_003442.2).
Other Names	ZNF177; PIGX; zinc finger protein 177
Accession No.	Swiss-Prot#:Q13360NCBI Gene ID:7730
Uniprot	Q13360
GeneID	100529215;7730;
Calculated MW	55kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse testis, using ZNF177 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only