

GASP2 Antibody

Catalog No: #45992

Package Size: #45992-1 50ul #45992-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GASP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	GASP2 Antibody detects endogenous levels of total GASP2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GASP2
Target Name	GASP2
Other Names	GPRASP2, GASP-2, GASP2
Accession No.	Swiss-Prot#: Q96D09NCBI Gene ID:100528062/114928
Uniprot	Q96D09
GeneID	100528062;114928;
Calculated MW	94kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis GASP2 using HeLa whole cell lysates

Note: This product is for in vitro research use only