

RALGPS1 Polyclonal Antibody

Catalog No: #28716



Package Size: #28716-1 50ul #28716-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RALGPS1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human RALGPS1 (NP_055451.1).
Other Names	RALGPS1; RALGEF2; RALGPS1A; RaI GEF with PH domain and SH3 binding motif 1
Accession No.	Swiss-Prot#:Q5JS13NCBI Gene ID:9649
Uniprot	Q5JS13
GeneID	9649;
Calculated MW	62kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RALGPS1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only