

TOB1 Polyclonal Antibody

Catalog No: #28721

Package Size: #28721-1 50ul #28721-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TOB1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TOB1 (NP_005740.1).
Other Names	TOB1; APRO5; APRO6; PIG49; TOB; TROB; TROB1; protein Tob1
Accession No.	Swiss-Prot#:P50616NCBI Gene ID:10140
Uniprot	P50616
GeneID	10140;
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TOB1 at 1:1000 dilution.

Background

This gene encodes a member of the transducer of erbB-2 /B-cell translocation gene protein family. Members of this family are anti-proliferative factors that have the potential to regulate cell growth. The encoded protein may function as a tumor suppressor. Alternate splicing results in multiple transcript variants.

Note: This product is for in vitro research use only