

DEXI Polyclonal Antibody

Catalog No: #28740



Package Size: #28740-1 50ul #28740-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DEXI Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human DEXI (NP_054734.2).
Other Names	DEXI; MYLE; Dexi homolog
Accession No.	Swiss-Prot#:O95424NCBI Gene ID:28955
Uniprot	O95424
GeneID	28955;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only