

SH3KBP1 binding protein 1 antibody

Catalog No: #23162

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SH3KBP1 binding protein 1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 162 and 447 of Human SHKBP1
Target Name	SH3KBP1 binding protein 1
Other Names	PP203
Accession No.	Swiss-Prot:Q8TBC3Gene ID:92799
Uniprot	Q8TBC3
GeneID	92799;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 76kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images




Sample(30 ug of whole cell lysate)

A: A431

B: H1299

7.5% SDS PAGE

Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded PC13 xenograft, using SHKBP1 antibody at 1: 100 dilution.

Note: This product is for in vitro research use only