

ZNF562 Polyclonal Antibody

Catalog No: #28750



Package Size: #28750-1 50ul #28750-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF562 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF562 (NP_001123503.1).
Other Names	ZNF562; zinc finger protein 562
Accession No.	Swiss-Prot#:Q6V9R5NCBI Gene ID:54811
Uniprot	Q6V9R5
GeneID	54811;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded human lung cancer using ZNF562 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human gastric cancer using ZNF562 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse liver using ZNF562 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using ZNF562 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only