

BATF3 Polyclonal Antibody

Catalog No: #28753

Package Size: #28753-1 50ul #28753-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BATF3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human BATF3 (NP_061134.1).
Other Names	BATF3; JDP1; JUNDM1; SNFT; basic leucine zipper ATF-like transcription factor 3
Accession No.	Swiss-Prot#:Q9NR55NCBI Gene ID:55509
Uniprot	Q9NR55
GeneID	55509;
Calculated MW	14kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:100

Images



Western blot analysis of extracts of rat spleen, using BATF3 at 1:1000 dilution.

Background

This gene encodes a member of the basic leucine zipper protein family. The encoded protein functions as a transcriptional repressor when heterodimerizing with JUN. The protein may play a role in repression of interleukin-2 and matrix metalloproteinase-1 transcription.

Note: This product is for in vitro research use only