

## VIPAS39 Polyclonal Antibody

Catalog No: #28757

Package Size: #28757-1 50ul #28757-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VIPAS39 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VIPAS39 (NP_071350.2).
Other Names	VIPAS39; C14orf133; SPE-39; SPE39; VIPAR; VPS16B; hSPE-39; spermatogenesis-defective protein 39 homolog
Accession No.	Swiss-Prot#:Q9H9C1NCBI Gene ID:63894
Uniprot	Q9H9C1
GeneID	63894;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

## Background

This gene encodes a protein involved in the sorting of lysosomal proteins. Mutations in this gene are associated with ARCS2 (arthrogryposis, renal dysfunction, and cholestasis-2). Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only