

# ZNF574 Polyclonal Antibody

Catalog No: #28760



Package Size: #28760-1 50ul #28760-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | ZNF574 Polyclonal Antibody                                        |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | IHC,IF                                                            |
| Species Reactivity    | Human,Mouse,Rat                                                   |
| Immunogen Description | Recombinant fusion protein of human ZNF574 (NP_073589.4).         |
| Other Names           | ZNF574; FP972; zinc finger protein 574                            |
| Accession No.         | Swiss-Prot#:Q6ZN55NCBI Gene ID:64763                              |
| Uniprot               | Q6ZN55                                                            |
| GeneID                | 64763;                                                            |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |


## Application Details

IHC 1:50 - 1:200IF 1:50 - 1:200


## Images




Immunofluorescence analysis of NIH-3T3 cells using ZNF574 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence analysis of U-2 OS cells using ZNF574 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunohistochemistry of paraffin-embedded human tonsil using ZNF574 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human lung using ZNF574 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human smooth muscle using ZNF574 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human breast cancer using ZNF574 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human stomach using ZNF574 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human uterus using ZNF574 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using ZNF574 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse pancreas using ZNF574 at dilution of 1:100 (40x lens).

## Background

Note: This product is for in vitro research use only