

LRFN3 Polyclonal Antibody

Catalog No: #28763



Package Size: #28763-1 50ul #28763-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LRFN3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LRFN3 (NP_078785.1).
Other Names	LRFN3; FIGLER1; SALM4; leucine rich repeat and fibronectin type III domain containing 3
Accession No.	Swiss-Prot#:Q9BTN0NCBI Gene ID:79414
Uniprot	Q9BTN0
GeneID	79414;
Calculated MW	66kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LRFN3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only