

SLC23A3 Polyclonal Antibody

Catalog No: #28783



Package Size: #28783-1 50ul #28783-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SLC23A3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human SLC23A3 (NP_653313.3).
Other Names	SLC23A3; E2BP3; SVCT3; Yspl1; solute carrier family 23 member 3
Accession No.	Swiss-Prot#:Q6PIS1NCBI Gene ID:151295
Uniprot	Q6PIS1
GeneID	151295;
Calculated MW	54kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse small intestine, using SLC23A3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only