

ZNF433 Polyclonal Antibody

Catalog No: #28787



Package Size: #28787-1 50ul #28787-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF433 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF433 (NP_001073880.1).
Other Names	ZNF433; zinc finger protein 433
Accession No.	Swiss-Prot#:Q8N7K0NCBI Gene ID:163059
Uniprot	Q8N7K0
GeneID	163059;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Immunofluorescence analysis of C6 cells using ZNF433 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.




Immunofluorescence analysis of NIH-3T3 cells using ZNF433 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.




Immunofluorescence analysis of U-2 OS cells using ZNF433 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.




Immunohistochemistry of paraffin-embedded rat brain using ZNF433 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human lung cancer using ZNF433 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human appendix using ZNF433 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human gastric cancer using ZNF433 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human uterine cancer using ZNF433 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse lung using ZNF433 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using ZNF433 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only