

ZNF549 Polyclonal Antibody

Catalog No: #28790



Package Size: #28790-1 50ul #28790-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF549 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ZNF549 (NP_001186224.1).
Other Names	ZNF549; zinc finger protein 549
Accession No.	Swiss-Prot#:Q6P9A3NCBI Gene ID:256051
Uniprot	Q6P9A3
GeneID	256051;
Calculated MW	68kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of Human placenta, using ZNF549 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only