

ZNF680 Polyclonal Antibody

Catalog No: #28791



Package Size: #28791-1 50ul #28791-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF680 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human ZNF680 (NP_848653.2).
Other Names	ZNF680; zinc finger protein 680
Accession No.	Swiss-Prot#:Q8NEM1NCBI Gene ID:340252
Uniprot	Q8NEM1
GeneID	340252;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IHC□1:50 - 1:200

Background

Note: This product is for in vitro research use only