

ASTL Polyclonal Antibody

Catalog No: #28793



Package Size: #28793-1 50ul #28793-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ASTL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ASTL (NP_001002036.3).
Other Names	ASTL; SAS1B; astacin-like metalloendopeptidase
Accession No.	Swiss-Prot#:Q6HA08NCBI Gene ID:431705
Uniprot	Q6HA08
GeneID	431705;
Calculated MW	46kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Raji cells, using ASTL at 1:1000 dilution.

Background

Note: This product is for in vitro research use only