

CHAD Polyclonal Antibody

Catalog No: #28794

Package Size: #28794-1 50ul #28794-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CHAD Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human CHAD (NP_001258.2).
Other Names	CHAD; SLRR4A; chondroadherin
Accession No.	Swiss-Prot#:O15335NCBI Gene ID:1101
Uniprot	O15335
GeneID	1101;
Calculated MW	40kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CHAD at 1:1000 dilution.

Background

Chondroadherin is a cartilage matrix protein thought to mediate adhesion of isolated chondrocytes. The protein contains 11 leucine-rich repeats flanked by cysteine-rich regions. The chondroadherin messenger RNA is present in chondrocytes at all ages.

Note: This product is for in vitro research use only