

PDHA2 Polyclonal Antibody

Catalog No: #28797



Package Size: #28797-1 50ul #28797-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PDHA2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse PDHA2 (NP_032837.1).
Other Names	PDHA2; PDHAL; pyruvate dehydrogenase (lipoamide) alpha 2
Accession No.	Swiss-Prot#:_NCBI Gene ID:18598
GeneID	:18598
Calculated MW	43kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using Pdha2 at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using PDHA2 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only